

Infectious Disease Collection Guide – Validated Sample Types



UTM Swab (MicroTest M6 Kit)

Micro-tipped and Traditional-Tipped swabs

- Respiratory Assays
 - COVID19, Flu A, Flu B, RSV A/B
 - Flu A, Flu B, RSV A/B
 - Adenovirus, Human
 - Metapneumovirus, Rhinovirus
 - Parainfluenza 1-4
- HSV 1+2/VZV



Aptima Multitest Swab Collection Kit

- STI Panel (females only)
 - *C. trachomatis*/*N. gonorrhoeae*
 - *T. vaginalis*
- *M. genitalium* (males and females)
- Women's Health (females only)
 - BV Assay
 - CV Assay



Aptima Urine Specimen Collection Kit (males and females)

- STI Panel
 - *C. trachomatis*/*N. gonorrhoeae*
 - *T. vaginalis*
 - *M. genitalium*

Note: urine must be transferred from a sterile container into the collection kit within 24 hours



SST

- HCV Viral Load
- HCV Genotyping



Orange Cap Fecal Swab

- Gastrointestinal Infection



Throat Swab (double swab)

- Strep Complete