

## Infectious Disease Collection Guide – Validated Sample Types



### Green/White Cap Swab

#### Green Cap- Thin Swab

Respiratory Assays (COVID 19, RPP targets, Strep Complete)

#### White Cap- Thick Swab

STI, *C. auris*, HSV 1+2/VZV

Note: White and green top tubes may be interchanged and tested on both sets of assays



### UTM Swab

- Respiratory Assays (COVID 19, RPP targets, Bordetella)
- Women's Health
- HSV 1+2/VZV



### Throat Swab (Double Swab)

- Strep Complete



### Urine with/without Boric Acid

- STI
- Urine Culture (Boric Acid)



### Orange Cap Fecal Swab

- Gastrointestinal Infection